

Business Reply  
Licence Number  
RRTA-YJHR-ZBEZ




dental  
solutions 

**Unit E35  
Champions Business Park  
Arrowe Brook Road  
UPTON, WIRRAL  
CH49 0AB**

Client:					
Address:					
Tel No:					
Patient:		Age:		M/F	
Indicate: <input type="checkbox"/> HEALTH SERVICE <input type="checkbox"/> INDEPENDENT <input type="checkbox"/> PRIVATE					
Indicate: <input type="checkbox"/> Crown(s) <input type="checkbox"/> Bridge(s) <input type="checkbox"/> Implant					
<input type="checkbox"/> Acrylic <input type="checkbox"/> Cobalt-Chrome <input type="checkbox"/> Orthodontic					
Bite:		Spec/ Tray			
Try-in		Re-try			
Shade/ Mould		Finish			
FOR LABORATORY USE ONLY					
Case Number:			Invoice Number:		
IN HOUSE INSPECTION					
Bite:		Spec/ Tray		Try-in	
Re-try		Metal Weight		Ceramic	
Approved for Manufacture by	Sign:	Date:	Approved for Release by	Sign:	Date:

### STAINS & CHARACTERISTICS



Surface Texture:  Heavy  
 Medium  
 Smooth

Glaze:  Low  
 Medium  
 High

**PRESCRIPTION FORM  
AND  
INSTRUCTION SHEET**

18	17	16	15	14	13	12	11	21	22	23	24	25	26	27	28
48	47	46	45	44	43	42	41	31	32	33	34	35	36	37	38

## NOTES Tel: 0151 678 9833

YOUR ATTENTION IS DRAWN TO THE FOLLOWING:  
 This is a custom made dental appliance that has been manufactured to satisfy the attributes, characteristics, properties and features specified by the client for the above patient. This dental appliance is intended for the exclusive use by this patient and conforms to the relevant essential requirements specified in Annex 1 of the Medical Devices Directive and the United Kingdom Medical Devices Regulations. THIS STATEMENT DOES NOT APPLY TO REPAIRS AND PPE MOUTHGUARDS.  
 STORING, HANDLING AND INSTRUCTIONS FOR USE:  
 It is recommended that before use this dental appliance is stored in a clean and safe environment that prevents it from coming into contact with materials, equipment, acids, alkalis or bleaches that could cause physical or chemical damage to the dental appliance. This dental appliance should not be subject to extremes of temperature during storage. Where applicable you should take care not to damage this dental appliance when removing it from the model. Where applicable, instructions on how to use or clean this dental appliance may be obtained from the prescriber.  
 THIS DENTAL APPLIANCE IS SUPPLIED IN AN UNSTERILISED STATE